## Learning about circulation... in a heartbeat!

#### Aim

Let's investigate what happens to our heart rate when we rest and exercise.

#### Method

Use the Curiscope 'Virtuali-Tee' app to take heart rate measurements and then note them down on the table.

Measure your heart rate every minute until it has returned to near the resting rate, this should be less than 10 minutes.

#### **Materials**

Timer or stop watch Curiscope 'Virtuali-Tee' app heart rate monitor


#### Results

| Time<br>(minutes) | Activity  | Heart rate (BPM) |
|-------------------|---|------------------|
| 0                 | Resting   |                  |
| 1                 | Complete<br>1 minute of exercise<br>(e.g. star jumps) |                  |
| 2                 |   |                  |
| 3                 |   |                  |
| 4                 |   |                  |
| 5                 |   |                  |
| 6                 |   |                  |
| 7                 | Recovering  |                  |
| 8                 |   |                  |
| 9                 |   |                  |
| 10                |   |                  |
| 11                |   |                  |
| 12                |   |                  |

www.curiscope.com

The ULTIMATE way to learn about the body





# The ULTIMATE way to learn about the body

www.curiscope.com

#### Discussion

- 1. What happened to your heart rate when you exercised?
- 2. Did your heart rate stay the same after exercising? Describe what happened.
- 3. Did you notice anything different about your breathing rate while you were doing the experiment? Describe what happened.

#### Conclusion

What did you discover about your heart rate when you rest, exercise and recover?

#### **Extension**

Can you devise an investigation to show the effects of different activities on heart rate? Does not need to be sports!

www.curiscope.com

The ULTIMATE way to learn about the body

## Group organisation – 4 students per group

**Student 1: The Human Subject** Does the exercise activity and has heart rate read as instructed. Student 2: The Timer Times the activity and informs other group members when to take measurements.

**Student 3: Heart Rate Measurer** Takes the heart rate of the Subject when instructed.

Warning: the heart rate monitor can get hot and the Subject must be free to move away from the light. **Student 4: The Recorder** Writes down measurements as taken by the Heart Rate Measurer.



www.curiscope.com

### Group organisation – 3 students per group

**Student 1: The Human Subject** Does the exercise activity and measures heart rate when instructed.


Warning: the heart rate monitor can get hot and the Subject must be free to move away from the light. Student 2: The Timer Times the activity and informs other group members when to take measurements.

**Student 3: The Recorder** Reads the heart rate and writes down

measurements as taken by the Heart Rate Measurer.

The ULTIMATE way to learn about the body

www.curiscope.com



-www.curiscope.com